
Electronic Supplementary Material (ESI) for RSC Advances. This journal is © The Royal Society of Chemistry 2015

Mwabcd230kD
150kD
80kD
60kDImage: Constraint of the second second

Figure S1 The result of SDS-PAGE analysis.



a, 10 uL of protein molecular weight marker; b, 10 uL of protein mixture solution containing BSA and Ova (each with a concentration of 1.0 mg/mL); c, 10 uL of mixture solution afteradsorption by the MIP; d, 10 uL of the mixture solution after adsorption by the NIP.