


**A ratiometric fluorescence sensor for HOCl based on FRET platform and application
in living cells**

Yan-Ru Zhang,^{[a]#} Zhi-Min Zhao,^{[b]#} Le Su,^[b] Jun-Ying Miao,^[b] and Bao-Xiang Zhao*^[a]

^aInstitute of Organic Chemistry, School of Chemistry and Chemical Engineering, Shandong University, Jinan 250100, P.R. China.

E-mail: bxzhao@sdu.edu.cn.

^bInstitute of Developmental Biology, School of Life Science, Shandong University, Jinan 250100, P.R. China.



Scheme S1 Synthesis of CRSH




Fig. S1 HRMS spectra of one fluorescent compound separated from the reaction of rhodamine thiohydrazide and HOCl.




Fig. S2 The red line is the donor (1 μM) and the green line is **CRSH** (1 μM) after addition of HOCl (4 μM) Condition: NaH_2PO_4 (0.05 M, pH = 5) : EtOH = 5 : 5 (v/v), λ_{ex} : 350 nm (slit widths: 10 nm/5 nm).

$$\text{Energy transfer efficiency (CRSH)} = 1 - F_{\text{DA}}/F_{\text{D}} = 1 - 46.3481/140.4140 = 67.0\%$$




Fig. S3 pH-dependent fluorescence intensity ratio changes of CRSH, Condition: [CRSH] = 1 μ M, NaH₂PO₄ (0.05 M, pH = 4-10) : EtOH (5 : 5, v/v), λ_{ex} : 350 nm (slit widths: 10 nm/10 nm), HOCl (4 μ M).




Fig. S4 (a) The ratio (I_{580}/I_{470}) of CRSH with addition of ROS/RNS. 1. free CRSH, 2. CRSH+HOCl, 3-10. CRSH+other ROS/RNS ($^1\text{O}_2$, H_2O_2 , $\text{HO}\cdot$, NO, NOOO^- , $^-\text{O}_2$, $t\text{-BuO}\cdot$, $t\text{-BuOOH}$)+HOCl. (b) Absorption spectra of CRSH toward HOCl (4 μ M), other ROS/RNS (20 μ M for $\text{HO}\cdot$, ONOO^- , NO, H_2O_2 , $t\text{-BuOOH}$, $t\text{-BuOO}\cdot$, $^1\text{O}_2$, $^-\text{O}_2$), Condition: Condition: [CRSH] = 1 μ M, NaH₂PO₄ (0.05 M, pH = 5) : EtOH (5 : 5, v/v).




Fig. S5 Absorption spectra of CRSH upon addition of HOCl (0–4 μM), Condition: [CRSH] = 1 μM , NaH_2PO_4 (0.05 M, pH = 5) : EtOH (5 : 5, v/v).




Fig. S6 Time-dependent fluorescence intensity ratio changes of CRSH, Condition: [CRSH] = 1 μM , NaH_2PO_4 (0.05 M, pH = 5) : EtOH (5 : 5, v/v), λ_{ex} : 350 nm (slit widths: 10 nm/10 nm).









Fig. S7 Viability of RAW264.7 cells after treatment with CRSH for 6 h (a), 12 h (b) at the different concentration (1, 5, 10 μM).







Sample Name: 1030-CRSH Position: P1-63 Instrument Name: Instrument 1 User Name:
 Inj Vol: -1 Inj/Position: Sample Type: Sample IRM Calibration Status: Success
 Data Filename: 1030-CRSH.d ACQ Method: 0103.m Comment: Acquired Time: 10/30/2015 2:58:49 PM

CRSH

