

1

Supplementary Information

2 **A novel colorimetric sensing platform for the detection of *S. aureus* with high sensitivity and**
3 **specificity**

4 Yun Zhang ^{a,1,*}, Shuyou Shi ^{a,1}, Jiajia Xing ^b, Wenqing Tan ^a, Chenguang Zhang ^a, Lin Zhang ^c, Huan
5 Yuan ^a, Miaomiao Zhang ^a, Jinjuan Qiao ^d

6

7 ^a Henan Key Laboratory of Immunology and Targeted Therapy, Henan Collaborative Innovation Center
8 of Molecular Diagnosis and Laboratory Medicine, School of Laboratory Medicine, Xinxiang Medical
9 University, Xinxiang 453003, PR China

10 ^b School of International Education, Xinxiang Medical University, Xinxiang 453003, PR China

11 ^c School of Innovation and Entrepreneurship, Xinxiang Medical University, Xinxiang 453003, PR China

12 ^d Department of Medical Laboratory, Weifang Medical University, Weifang 261053, PR China

13


14 ¹ these authors contributed equally to this work

15

16 *Corresponding authors: Yun Zhang, Tel: +86 373 3029977; Fax: + 86 373 3029977; Email:

17 zhangyun0126@126.com

18




19

20 **Fig. S1** (A) The image of blue bacteria pellet of Cy5.5 binding *S. aureus* after centrifuge. (B)

21 Fluorescence micrographs of Cy5.5-stained *S. aureus*. Bar size: 50 μm

22



23

24 **Fig. S2** Effects of experimental conditions on signals of the system. Three independent measurements

25 were taken from three individual preparations for each condition. Error bars indicated the standard

26 deviations


27

Table S1 Standard deviation test data of negative control

No.	Absorbance at 450 nm (a.u.)
1	0.049
2	0.05
3	0.051
4	0.051
5	0.052
6	0.052
7	0.052
8	0.052
9	0.052
10	0.057

Table S2 Relative standard deviation test data for reproducibility


No.	Absorbance at 450 nm (a.u.)
1	1.218
2	1.22
3	1.222
4	1.241
5	1.241
6	1.267
7	1.267
8	1.293
9	1.334
10	1.250
11	1.356



34

35 **Fig. S3** Effect of pH on the performance of the proposed method. Red bars indicated absorbance intensity
36 (450 nm) (from left to right) for 2.0×10^5 CFU/mL *S. aureus* in buffer solutions over a range of pH values
37 3-9. Error bars indicated the standard deviations


38



39

40 **Fig. S4** Effect of salt concentration of buffer solution on the performance of the proposed method. Red
41 bars indicated absorbance intensity (450 nm) (from left to right) for 2.0×10^5 CFU/mL *S. aureus* in buffer
42 solutions over a range of PBS concentrations (0.01 – 100 mmol/L). Error bars indicated the standard
43 deviations


44



45

46 **Fig. S5** Calibration curve of absorbance intensity (450 nm) along with *S. aureus* concentration under the
47 optimal conditions in orange juice (the sample was 5-time diluted). Three independent measurements
48 were taken from three individual preparations for each condition. Error bars indicated the standard
49 deviations

50




51

52 **Fig. S6** Calibration curve of absorbance intensity (450 nm) along with *S. aureus* concentration under the

53 optimal conditions in spring water. Three independent measurements were taken from three individual

54 preparations for each condition. Error bars indicated the standard deviations

55



56

57 **Fig. S7** Calibration curve of absorbance intensity (450 nm) along with *S. aureus* concentration under the
58 optimal conditions in human urine. Three independent measurements were taken from three individual
59 preparations for each condition. Error bars indicated the standard deviations

60


61 **Table S3** Analytical results of *S. aureus* in real samples obtained using the proposed method and by

62 ELISA (n=3)

Samples	Spiked (CFU mL⁻¹)	Found (CFU mL⁻¹)	RSD (%)	ELISA (CFU mL⁻¹)	RSD (%)
Orange juice ^a	1.3 × 10 ⁴	1.4 × 10 ⁴	6.32	1.2 × 10 ⁴	2.31
	5.0 × 10 ⁴	5.5 × 10 ⁴	7.21	5.4 × 10 ⁴	4.42
Spring water	1.3 × 10 ⁴	1.2 × 10 ⁴	6.87	1.1 × 10 ⁴	2.21
	5.0 × 10 ⁴	5.2 × 10 ⁴	4.35	4.7 × 10 ⁴	2.87
Human urine ^b	1.3 × 10 ⁴	1.2 × 10 ⁴	3.52	1.3 × 10 ⁴	2.39
	5.0 × 10 ⁴	5.0 × 10 ⁴	1.75	5.5 × 10 ⁴	1.09

63 ^a The sample was 5 times diluted. ^b The sample was 10 times diluted.

64



65

66 **Fig. S8** Different strains of *S. aureus* detected by the proposed method. Red bars indicated absorbance

67 intensity (450 nm) (from left to right) for 2.0×10^5 CFU/mL *S. aureus* (ATCC 139843), *S. aureus* (ATCC

68 25923) and *S. aureus* (ATCC 12598). Error bars indicated the standard deviations