

Supplementary

Figure legend

Fig. S1. (a) FTIR of pure TGA, CdTe QD_(D) and CdTe/CdS QD_(A), (b) FT-IR spectra of CdTe QD_(D), QD_(D)* conjugated with mAb and CdTe/CdS QD_(A), QD_(A)* conjugated with pAb

Fig. S2- Dependence of $(ahn)^2$ of CdTe core and CdTe/CdS coreshell thin films upon the incident photon energy ($h\nu$).

Fig. S3. (a) Absorption spectra and fluorescence of QD(D) alone and QD(D)-mAb, QD(A) alone with QD(A)-pab. (b) Fluorescence spectra of QD_(D) and (b) QD_(D)-mAb, QD_(A) and QD_(A)-pAb.

1 **Fig. S4.** DLS measurement of (a), QD(D) alone (b)QD(D)-mAb, (c)QD(A) alone and (d) QD(A)-pab

Fig. S5- Stability of the FRET based method for detection of NMP22.

Fig. S6- Study the effect of interferences on FRET based method for detection NMP22.

Table legend

Table S1- The comparison with other methods for NMP 22.

Table S2- determination of NMP 22 in urine real samples.





Fig. S1. (a) FTIR of pure TGA, CdTe QD_(D) and CdTe/CdS QD_(A), (b) FT-IR spectra of CdTe QD_(D), QD_(D)* conjugated with mAb and CdTe/CdS QD_(A), QD_(A)* conjugated with pAb



1 Fig. S2. Dependence of $(ahn)^2$ of CdTe core and CdTe/CdS coreshell thin films upon the incident photon energy ($h\nu$)




Fig. S3. (a) Absorption spectra and fluorescence of QD(D) alone and QD(D)-mAb, QD(A) alone with QD(A)-pAb. (b) Fluorescence spectra of QD(D) and (b) QD(D)-mAb, QD(A) and QD(A)-pAb





Fig. S6. Study the effect of interferences on FRET based method for detection NMP22

Table S1. The comparison with other methods for NMP 22

Detection technique	Materials used	Linear range	LOD pg.mL^{-1}	References
Sandwich electrochemical immunoassay		(1.2–200) ng.mL^{-1}	500	³¹
Alere NMP22 Test reagent kits	HRP-Ab-NMP 22/Fe ₃ O ₄ /Au/CoPc CME	(4–232) ng.mL^{-1}	-	³²
ECL immunoassay	-	(0.05–2) ng.mL^{-1}	10	¹
FRET immunosensor	g-C ₃ N ₄ @Au-Ab-NMP22 QD _(D) -mAb-NMP22	(2–22) pg.mL^{-1}	0.05	This work

Table S2. determination of NMP 22 in urine real samples

#	FRET method (pg.mL^{-1})	Standard method (pg.mL^{-1})	Recovery%	RSD% n=5
1	18.87	18	104.84	2.35
2	12.01	12	100.07	1.85
3	6.30	6	105.00	3.00
4	3.03	3	101.00	2.56

