

## Supplementary Materials

### Contents

1. Reference substances for fatty acid composition analysis .....	1
2. Identification of volatile components of SSO .....	3
2.1 Methods.....	3
2.2 Results.....	3
3. Western blots uncropped images.....	5
3.1 GLUT9.....	5
3.2 URAT1 .....	6
3.3 OAT1.....	7
3.4 $\beta$ -actin for GLUT9, URAT1 and OAT1 .....	8
3.5 Cytosol Keap1.....	9
3.6 Cytosol Nrf2.....	10
3.7 $\beta$ -actin for cytosol Keap1 and Nrf2.....	11
3.8 Nuclear Nrf2 .....	12
3.9 Histone H3 .....	13
3.10 TXNIP.....	14
3.11 NLRP3.....	15
3.12 ASC .....	16
3.13 Caspase1 .....	17
3.14 $\beta$ -actin for TXNIP, NLRP3, ASC and Caspase1 .....	18

## 1. Reference substances for fatty acid composition analysis

Reference substances used in the fatty acid composition analysis of *Sonneratia apetala* seed oil (SSO) were shown in Fig.1 and Tab.1. GLC MIXTURE 617 (C4 - C24, purity  $\geq 98.0\%$ ) were purchased from Nu-Chek Prep, Inc (Minnesota, USA).




Fig.1 Total ion chromatogram of reference substances for fatty acid composition analysis

Tab.1 Reference substances for fatty acid composition analysis of *Sonneratia apetala* seed oil

Number	Retention time (min)	Molecular weight	Identification	Wt. %
1	5.879	88.11	Butyrate, methyl ester	2
2	6.684	116.16	Hexanoate, methyl ester	2
3	8.024	144.21	Octanoate, methyl ester	2
4	9.922	172.26	Decanoate, methyl ester	4
5	11.065	186.29	Hendecanoate, methyl ester	2
6	12.357	200.32	Dodecanoate, methyl ester	4
7	13.809	214.34	Tridecanate, methyl ester	2
8	15.447	228.37	Tetradecanoate, methyl ester	2
9	16.766	226.35	(Z)-9-Tetradecenoate, methyl ester	2
10	17.281	242.40	Pentadecanoate, methyl ester	2
11	18.754	240.38	(Z)-10-Pentadecenoate, methyl ester	2
12	19.315	256.42	Hexadecanoate, methyl ester	4
13	20.616	254.41	(Z)-9-Hexadecenoate, methyl ester	2
14	21.508	270.45	Heptadecanoate, methyl ester	4
15	22.900	268.43	(Z)-10-Heptadecenoate, methyl ester	2

16	23.843	284.48	Octadecanoate, methyl ester	4
17	24.677	282.46	(E)-9-Octadecenoate, methyl ester	2
18	25.104	282.46	(Z)-9-Octadecenoate, methyl ester	4
19	26.155	280.44	(E,E)-9,12-Octadecadienoate, methyl ester	2
20	27.120	280.44	(Z,Z)-9,12-Octadecadienoate, methyl ester	2
21	28.598	278.44	(Z,Z,Z)-6,9,12-Octadecatrienoate, methyl ester	4
22	28.752	312.53	Eicosanoate, methyl ester	4
23	29.549	278.44	(Z,Z,Z)-9,12,15-Octadecatrienoate, methyl ester	2
24	30.044	310.52	(Z)-11-Eicosenoate, methyl ester	2
25	31.237	326.56	Heneicosanoate, methyl ester	2
26	32.166	308.53	(Z,Z)-11,14-Eicosadienoate, methyl ester	2
27	33.661	306.48	(Z,Z,Z)-8,11,14-Eicostrienoate, methyl ester	2
28	33.785	340.58	Docosanoate, methyl ester	4
29	34.640	306.48	(Z,Z,Z)-11,14,17-Eicostrienoate, methyl ester	2
30	34.754	304.47	(Z,Z,Z,Z)-5,8,11,14-Eicosatetraenoate, methyl ester	2
31	35.084	338.57	(Z)-13-Docosenoate, methyl ester	2
32	36.240	336.55	(Z,Z)-13,16-Docosadienoate, methyl ester	2
33	37.206	354.61	Tricosanoate, methyl ester	2
34	37.306	302.45	(Z,Z,Z,Z,Z)-5,8,11,14,17-Eicosapentaenoate, methyl ester	2
35	38.729	368.64	Tetracosanoate, methyl ester	4
36	39.961	332.51	(Z,Z,Z,Z)-7,10,13,16-Docosatetraenoate, methyl ester	2
37	40.026	366.62	(Z)-15-Tetracosenoate, methyl ester	2
38	40.970	330.50	(Z,Z,Z,Z,Z)-4,7,10,13,16-Docosapentaenoate, methyl ester	2
39	42.550	330.50	(Z,Z,Z,Z,Z)-7,10,13,16,19-Docosapentaenoate, methyl ester	2
40	43.573	328.49	(Z,Z,Z,Z,Z,Z)-4,7,10,13,16,19-Docosahexaenoate, methyl ester	2

## 2. Identification of volatile components of SSO

### 2.1 Methods

SSO sample (100  $\mu\text{L}$ ) was extracted with 300  $\mu\text{L}$  of n-hexane and 1 mL of acetonitrile under ultrasound for 20 min. Following centrifugation at  $14000 \times g$  for 10 min ( $4^\circ\text{C}$ ), supernatant was obtained, and the remains was extracted with 1 mL of acetonitrile for another 2 min. The upper extract was also harvested after centrifugation for 10 min ( $14000 \times g$ ,  $4^\circ\text{C}$ ) and mixed with the previously collected supernatant for the GC-MS profiling. The yield of volatile components of SSO was 9.56%.

A total of 5  $\mu\text{L}$  sample was separated with DB-wax GC-Column (30 m  $\times$  0.25 mm, 0.25  $\mu\text{m}$ ), carried by mixed gas of nitrogen and helium at flow rate of 1 mL/min. The temperature of column injector was maintained at  $250^\circ\text{C}$  and split ratio was set at 30:1. Oven temperature was held at  $50^\circ\text{C}$  for 2 min and elevated to  $280^\circ\text{C}$  at rate of  $10^\circ\text{C}/\text{min}$  and held for 20 min. The mass scan parameters were set as followings: electron impact ionization voltage 70 eV; ion source temperature  $250^\circ\text{C}$ ; mass range 30 ~ 450 m/z.

### 2.2 Results

As shown in Fig.2 and Tab.2, the dominant volatile constituents of SSO were free fatty acids and related derivatives including linoleic acid, linoleic acid ethyl ester, (E)-2-heptenal, (E, E)-2,4-decadienal, palmitic acid, oleic acid, etc. Additionally, small amounts of terpenes including squalene, limonene and  $\alpha$ -phellandrene were also found in SSO.




Fig.2 Total ion chromatogram of volatile organic components analysis




Tab. 2 Volatile organic components of *Sonneratia apetala* seed oil

Number	Retention time (min)	Molecular weight	Molecular formula	Indetification	% composition
1	6.829	136.23	C <sub>10</sub> H <sub>16</sub>	Limonene	0.12
2	6.982	136.23	C <sub>10</sub> H <sub>16</sub>	$\alpha$ -Phellandrene	0.12
3	8.734	112.17	C <sub>7</sub> H <sub>12</sub> O	(E)-2-Heptenal	5.25
4	14.858	152.23	C <sub>10</sub> H <sub>16</sub> O	(E, E)-2,4-Decadienal	4.74
5	19.148	284.48	C <sub>18</sub> H <sub>36</sub> O <sub>2</sub>	Hexadecanoic acid ethyl ester	2.47
6	21.121	310.51	C <sub>20</sub> H <sub>38</sub> O <sub>2</sub>	(E)-9-Octadecenoic acid ethyl ester	1.23
7	21.541	308.50	C <sub>20</sub> H <sub>36</sub> O <sub>2</sub>	Octadecadienoic acid ethyl ester	9.25
8	24.519	256.42	C <sub>16</sub> H <sub>32</sub> O <sub>2</sub>	n-Hexadecanoic acid	4.62
9	25.513	410.70	C <sub>30</sub> H <sub>50</sub>	Squalene	0.95
10	25.973	284.48	C <sub>18</sub> H <sub>36</sub> O <sub>2</sub>	Octadecanoic acid	2.25
11	26.207	282.46	C <sub>18</sub> H <sub>34</sub> O <sub>2</sub>	Octadecenoic acid	2.89
12	26.614	280.45	C <sub>18</sub> H <sub>32</sub> O <sub>2</sub>	(Z, Z)-9,12-Octadecadienoic acid	12.64

### 3. Western blots uncropped images




#### 3.1 GLUT9 (The band in the red frame is the target band)

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO






**3.2 URAT1 (The band in the red frame is the target band)**

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO



### 3.3 OAT1 (The band in the red frame is the target band)




Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO








**3.4  $\beta$ -actin for GLUT9, URAT1 and OAT1 (The band in the red frame is the target band)**

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO






### 3.5 Cytosol Keap1 (The band in the red frame is the target band)

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO






### 3.6 Cytosol Nrf2 (The band in the red frame is the target band)

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO






**3.7  $\beta$ -actin for cytosol Keap1 and Nrf2 (The band in the red frame is the target band)**

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO






**3.8 Nuclear Nrf2 (The band in the red frame is the target band)**

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO






### 3.9 Histone H3 (The band in the red frame is the target band)

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO






**3.10 TXNIP (The band in the red frame is the target band)**

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO



### 3.11 NLRP3 (The band in the red frame is the target band)




Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO








**3.12 ASC (The band in the red frame is the target band)**

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO



### 3.13 Caspase1 (The band in the red frame is the target band)

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO



**3.14  $\beta$ -actin for TXNIP, NLRP3, ASC and Caspase1(The band in the red frame is the target band)**

Intact, Vehicle, BZM, FBX, 50 mg/kg of SSO, 100 mg/kg of SSO, 200 mg/kg of SSO

