

Supporting Information

Microdroplet Extraction Assisted Ultrasensitive Raman Detection in Complex Oil

Lingjun Zhang,^{‡a} Jianfeng Sun,^{‡a} Hong He,^a Yingzhou Huang,^{*a} Haiyang Shi^{*c} and Weigen Chen^b

^a State Key Laboratory of Coal Mine Disaster Dynamics and Control, Chongqing Key Laboratory of Soft Condensed Matter Physics and Smart Materials, College of Physics, Chongqing University, Chongqing, 400044, China. E-mail: yzhuang@cqu.edu.cn (Y. Huang)

^b State Key Laboratory of Power Transmission Equipment & System Security and New Technology, Chongqing University, Chongqing 400044, China.

^c School of Mechatronics and Vehicle Engineering, Chongqing Jiaotong University Chongqing 400074, China. E-mail: shihaiyang@cqu.edu.cn (H. Shi).

† Electronic supplementary information (ESI) available. See DOI: xxx

‡ These authors contributed equally to this work.

List of Contents

Fig. S1 (a) The geometry of the microfluidic channel. The serpentine channel and the microdroplet Raman detection channel were 150 and 800 μm wide, respectively. The channel height was 100 μm . (b) Schematic of experimental setup of fluorescence measurements, comprising the microfluidic chip, fluorescence microscope and two syringe pumps.

Fig. S2 A schematic of a home-built Raman spectrometer used for microdroplet interrogation.




Fig. S1 (a) The geometry of the microfluidic channel. The serpentine channel and the microdroplet Raman detection channel were 150 and 800 μm wide, respectively. The channel height was 100 μm . (b) Schematic of experimental setup of fluorescence measurements, comprising the microfluidic chip, fluorescence microscope and two syringe pumps.




Fig. S2 A schematic of a home-built Raman spectrometer used for microdroplet interrogation.

Table S1. Sensitivity analysis between HPLC and microdroplet Raman.

Concentration of furfural in samples (ppm)	HPLC		Microdroplet extraction assisted real-time Raman spectroscopy microfluidic chip				
	Detection result (ppm)	Deviation (%)	Detection result (ppm)			Deviation (%)	Coefficient of standard deviation (%)
			1	2	3		
1.00	1.07	7.00	1.01	0.95	1.07	4.72	5.15
2.00	1.92	4.00	1.77	2.20	1.96	7.80	8.83
3.30	3.48	5.45	3.32	3.07	3.59	5.37	6.34