

Supporting Information

Cholesteroled polymer (Chol-b-Lys)-based nanoparticles (CL-NPs) confer antibacterial efficacy without resistance

Huan-Huan Ding,^a Jia-Zhu Chigan,^a Jian-Bin Zhen,^a Lu Liu,^a Yin-Sui Xu,^a Cheng Chen,^a and Ke-Wu Yang^{*a}

Proton Nuclear Magnetic Resonance (¹H NMR)

The NMR spectrum was recorded on a Bruker AV 400 MHz spectrometer, using tetramethylsilane as an internal standard and DMSO-d₆ as solvent. The sample concentration was ca. 10 mg/mL.




Fig. S1 ¹H NMR spectrum of L-Lys-Cbz-NCA in DMSO-d₆




Fig. S2 ^1H NMR spectrum of Chol- NH_2 in $\text{CCl}_3\text{-d}$




Fig. S3 ^1H NMR spectrum of Chol-*b*-Lys-Cbz in DMSO-d_6




Fig. S4 ^1H NMR spectrum of Chol-*b*-Lys in DMSO-d_6

Dynamic light scattering apparatus (DLS)

Measurement of particle size Distribution Polymer Micellar solution was carried out by dynamic light scattering apparatus (DLS) on a Malvern Zetasizer Nano-ZS at a fixed scattering angle of 90° . Measured with polymer concentration of 1 mg/mL for sample preparation.




Fig. S5 DLS studies of CL-NPs (by intensity) in water

Formulation of M9 medium

MgSO₄ (1 M), CaCl₂ (1 M) and ZnSO₄ (1 M) were dissolved in double distilled water (10 mL), respectively, and autoclaved for later use. 5×M9 salt solution, including Na₂PO₄·7H₂O, KH₂PO₄, NaCl and NH₄Cl, was dissolved in 200 ml of double distilled water and sterilized at 121 °C for 15 min. Then a 20% glucose solution was prepared with sterilized at 115 °C for 15 min.

Aseptic preparation of M9 medium (1000 mL): 5×M9 salt solution (200 mL), 1 M MgSO₄ (2 mL), 1 M CaCl₂ (0.1 mL), 1 M ZnSO₄ (0.1 mL) and 20% glucose solution (20 mL) were mixed, and the sterilized double distilled water was added to the mixture to reach a total volume of 1000 mL.

Scanning electron microscopy (SEM)

The morphology of the microorganisms before and after treatment with polymer was observed using a field emission SEM (SU8010) operated at an accelerating voltage of 3.0 keV. Samples were treated with gold before observation.

Zeta potential

Zeta potentials were measured on a Water Nano-ZS 90 Nanosizer (Malvern Instrument) at a fixed scattering angle of 90° at room temperature.