

1 Porous Organic Polymer Nanosphere-Based Fluorescent Biosensing Platform for
2 Simultaneous Detection of Multiplexed DNA via Electrostatic Attraction and π - π
3 Stacking Interactions

4 *Yujie Sun*[#], *Zhenzhong Lu*[#], *Wenlin Ma*, *Rui Wang*, *Chengwu Zhang*, *Jinhua Liu*^{*}

5 ,

6

7

8

9

10

11

12

13

14

15

16

17


18

19

20

21

22



23

24 **Fig. S1.** (a) The picture of POP powder; (b) TGA of POP.

25

26

27

28

29

30

31

32

33

34

35


36

37

38

39

40



41

42 **Fig. S2.** (a) Fluorescence emission spectra of P1, P1+POP and P1+ (POP/Mg²⁺); (b)

43 The zeta potentials of POP+H₂O and POP/Mg²⁺+H₂O; (C) Normalization intensity of

44 P1, P1+POP, P1+(POP/Metal ion) for Mg²⁺, Cd²⁺, Fe²⁺ and Fe³⁺, respectively.

45

46

47

48

49

50

51

52

53

54


55

56

57

58

59



60

61 **Fig. S3.** (a) Fluorescence emission spectra and quenching efficiency (b) of P1 (25 nM)

62 with different concentrations of POP (0-50 mg/mL).

63

64

65

66

67

68

69


70

71

72

73

74



75 **Fig. S4.** Fluorescence emission spectra of P1 (a), P2 (b), and corresponding quenching
 76 efficiency (c) with different concentrations of POP (0-60 mg/mL).

77

78

79

80

81

82

83

84

85

86

87


88

89

90

91

92



93 **Fig. S5.** Fluorescence intensity of P1 (a), P2 (b) and P3 (c) change with time in the
94 presence of POP and Target DNA successively.

95

96

97

98

99

100

101

102

103

104

105

106


107

108

109

110

111



112

113 **Fig. S6.** Fluorescence intensity of P1 change with time upon the introduction of
114 different quenchers (GO and POP, 30 $\mu\text{g/mL}$) and T1 (200 nM).

115

116

117

118

119

120

121

122

123


124

125

126

127

128



129

130 **Fig. S7.** (a) Fluorescence spectra of probes (P1, P2, and P3) under different conditions

131 (P, P+POP, P+POP+T), respectively. (b) Quenching efficiency of P (P1, P2 and P3)

132 caused by POP with different concentrations (0, 10, 20, 25, 30, 50, 70 and 90 $\mu\text{g/mL}$).

133

134

135

136

137

138

139


140

141

142

143

144



145

146 **Fig. S8.** Fluorescence emission spectra of P1 (a), P2 (b) and P3(c) change with different

147 concentrations of POP (0-90 mg/mL).

148

149

150

151

152

153

154

155

156

157

158

159


160

161

162

163

164



165

166 **Fig. S9.** Fluorescence spectra of probes (P1, P2 and P3) in one system toward T1 (a),
 167 T2 (b) and T3 (c); (d) Specificity of probes (P1, P2, P3) for targets (T1, T2, T3). B (1-
 168 3), blank; R (1-3), random sequence for P (1-3); M (1-3), base mismatched sequence
 169 for P (1-3), respectively.

170

171

172

173

174

175

176

177

178

179 **Table S1.** DNA sequences used in this work.

Oligonucleotide	Sequence (5-3)*
P1	CAGTTACATTCTCCCAGTTGATT-FAM
P2	ACCTGGGGGAGTATTGCGGAGGAAGGT-Cy3
P3	AGACTCTTGAGTTCTCAGTATG-Texas Red
T1	AATCAACTGGGAGAATGTA ACTG
T2	ACCTTCCTCCGCAATACTCCCCCAGGT
T3	CATACTGAGAACTCAAGAGTCT
M1	AATCAACTG <u>T</u> GAGAATGTA ACTG
M2	ACCTTCCTCCGCA <u>T</u> TACTCCCCCAGGT
M3	CAT <u>A</u> TTGAGAACTCAAGAGTCT
R	ACCTGGGGGAGTATTGCGGAGGAAGGT

180

181 * The mutation base is indicated by underline

182

183

184

185

186

187

188

189

190

191

192

193 **Table S2.** Comparison of different fluorescent DNA sensors.

Materials	Sensitivity	Detection time	Fluorescent reporter	Multiplexed detection	Reference
SWNT	4 nM	several hours	FAM	Not reported	[39]
AuNP	nM	minutes	Rhodamine 6G	Yes	[40]
Graphene oxide (GO)	100 pM	≈1 min	FAM, Cy3, ROX	Yes	[41]
GO (premixing)	≈10 nM	30 min	FAM	Not reported	[42]
GO	2 nM	30 min	FAM	Not reported	[43]
GO	0.5 nM	60-90 min	Silver	Yes	[44]
GO	75 nM	30 min	Graphene quantum dots	Not reported	[45]
g-C ₃ N ₄	81 pM	Not reported	FAM, ROX	No	[46]
MOFs(UiO-66-NH ₂)	10 nM	20 min	FAM	Not reported	[47]
MOFs(Cu(H ₂ dtoa))	3 nM	4 h	FAM	Not reported	[48]
POP	50 pM	25 min	FAM, Cy3, Texas Red	Yes	This work

194

195

196

197

198

199

200

201

202

203

204

205

206

207

208

209 **Table S3.** The quenching efficiency and recovery of different detected system.

System	Target DNA	QE (%)	Recovery (%)
Single chain	T1	94.43	38.37
	T2	83.12	55.90
	T3	85.80	53.04
Double chains	T1&T2	87.61	28.70
		68.18	83.09
	T2&T3	70.85	78.29
		48.92	70.34
	T1&T3	93.25	24.55
	86.92	29.68	
Three chains	T1&T2&T3	92.75	31.18
		87.62	28.64
		89.63	25.05

210

211

212