Electronic Supplementary Material (ESI) for Analyst. This journal is © The Royal Society of Chemistry 2022

Supporting Information

Hairpin-functionalized DNA tetrahedral for miRNA imaging in living cells via self-assembly to form dendrimers

Wenjing Zhao‡, Yao Jiang‡*, Huimin Zhou, Shusheng Zhang*

Shandong Provincial Key Laboratory of Detection Technology for Tumor Markers, Collaborative Innovation Center of Tumor Marker Detection Technology, Equipment and Diagnosis-Therapy Integration in Universities of Shandong, College of Chemistry and Chemical Engineering, Linyi University, Linyi 276005, P. R. China

Table S1. Sequences of Oligonucleotides Used in This Study

Name	Sequence (5'-3')	
H1	GACGACTAATAAGATTCAACATCAG/iBHQ1-dT/CTGATAAGCGCGCTCTCATAGGCGCTTATCAGACT-FAM	
H2	TGTCCTAATTAGAATATAAGCGCTCCGAGCCGGTCGAAATGAGAGCGCGCTTATCAGACACCGGCTCGGAG	
A1	ATCTTATTAGTCGTCTTCTCCAAGTCTACCTTGCTACACGACGCCATAGTAAAAAAAA	
A2	ATTCTAATTAGGACATTCACAGCAGTGCAGTGTAGCAAGAGGCGAGGGTCC	
A3	ATCTTATTAGTCGTCTTCACTGCTGTGAAACACTACGTGTCTACTATGGCG	
A4	ATTCTAATTAGGACATTAGACTTGGAGGCCACGTAGTGTTTTGGACCCTCGC	
miRNA-21	UAGCUUAUCAGACUGAUGUUGA	
miRNA-141	UAACACUGUCUGGUAAAGAUGG	
Let-7a	UGAGGUAGGUUGUAUAGUU	
Let-7b	UGAGGUAGGUUGUGUGUU	
Let-7c	UGAGGUAGGUUGUAUGGUU	
MUC1 aptamer	GGCTATAGCACATGGGTAAAACGACTTTTTTTTTT	





Figure S1. In vitro fluorescence response of free-H1. Schematic and fluorescence emission spectra of free-H1 with various concentrations of miRNA-21.



o fluorescence response of free-CHA. Schematic and fluorescence emission spectra of free-CHA with various concentrations of miRNA-21.




Figure S3. (a) RT-PCR measurements of miRNA-21 levels in different concentrations. (b) Measurements of miRNA-21 levels in MCF-7 cells by the HDFT and RT-PCR.




Figure S4. PPAGE characterization of (a) HFDT with different concentrations of miRNA; the stability of the (b) HFDT and (c) target-triggered assembly dendrimers structure.




Figure S5. Confocal fluorescence images of MCF-7 cells after incubation with HFDT at different times.




Figure S6. Flow cytometric quantification of (a) L02 cells and (b) MCF-7 cells treated with PBS and HDFT in different times.

Table S2. The summary of recent miRNA-21 analytical methods in comparison with our strategy.

Analysis method	Detection Technique	LOD	Reference
DNA tetrahedrons based fluorescent biosensor	Fluorescence	620 pM	1
Strand displacement reaction	Photoluminescence	4.5 pM	2
Pyrenelabelled Y-shape DNA nanoprobe based Fluorescent biosensor	Fluorescence	200 pM	3
Accelerated-DzFN	Fluorescence	1.2 pM	4
Hairpin-functionalized DNA tetrahedron (HFDT)	Fluorescence	0.8 pM	This work

References

- 1. W. Zhou, D. Li, C. Xiong, R. Yuan and Y. Xiang, ACS Appl. Mater. Inter., 2016, 8, 13303-13308.
- 2. X. He, T. Zeng, Z. Li, G. Wang and N. Ma, Angew. Chem. Int. Ed. Engl., 2016, 55, 3073-3076.
- 3. Z. Chen, J. Lu, F. Xiao, Y. Huang, X. Zhang and L. Tian, Chem. Commun., 2020, 56, 1501-1504.
- 4. Y. Wu, H. Meng, J. Chen, K. Jiang, R. Yang, Y. Li, K. Zhang, L. Qu, X. Zhang and Z. Li, Chem. Commun., 2020, 56, 470-473.