

Supporting information for

**Magnet-actuated microfluidic array chip for high-throughput pretreatment,
amplification and detection of multiple pathogens**

Mengfan Zhou^{§, a}, Huiying Su^{§, b, a}, Bangfeng Wang^a, Chao Wan^a, Wei Du^a,
Peng Chen^{*, a}, Xiaojun Feng^{*, a}, Bi-Feng Liu^a

^aThe Key Laboratory for Biomedical Photonics of MOE at Wuhan National Laboratory for Optoelectronics - Hubei Bioinformatics & Molecular Imaging Key Laboratory, Systems Biology Theme, Department of Biomedical Engineering, College of Life Science and Technology, Huazhong University of Science and Technology, Wuhan 430074, China

^bSchool of Biological engineering, Huainan Normal University, Huainan, Anhui 232038, China

§These authors contributed equally to this work.

* Corresponding author

Add: Systems Biology Theme, Department of Biomedical Engineering

College of Life Science and Technology

Huazhong University of Science and Technology, Wuhan 430074, China

Email: xfeng@mail.hust.edu.cn

gwchenpeng@mail.hust.edu.cn

Tel: +86-27-87793180

Fax: +86-27-87792170

Table S1 Sequences of the specific genes used for pathogen detection in this paper

Bacteria	Target Gene	GeneBnak Accession	Sequences (5'-3')
<i>Escherichia coli</i>	<i>mal B</i>	J01648	GCCATCTCCTGATGACGCATAGTCAGC
			CCATCATGAATGTTGCTGTCGATGACA
			GGTTGTTACAAAGGGAGAAGGGCATG
			GCGAGCGTACAGCTGCAAAATGTAAC
			GAAAGCCTGGGGCGAGGTCGTGGTAT
			CGAAAGATATCAATCTCGATATCCAT
			GAAGGTGAATTCGTGGTGTTTGTCGG
			ACCGTCTGGCTGCGGTAAAT
			GGCGATATTGGTGTTTATGGGGTCGTT
			CTACATTGACAGAATCCTCAGTTTTTC
			AACGTTTCCTGCGGTACTGTTAATTAC
			CACGCTCTTTCGTCTGGCATTATCGAT
			CAGTACCAGTCGTCTTATCTTGATTGA
			AGCCGATGCCGGTGAAATTATCGCCA
			CGTTCGGGCAATTCGTTATTGGCGATA
GCCTGGCGGTGGGTTTTGTTGTCTTCT			
CTATTGTCACCGTGGTCCAGTTTATCG			
			TT




Figure S1. The Standard curve of bacteria. (a) *Salmonella typhimurium*. (b) *Escherichia coli*.





Figure S2. (a) Miniature peltier heater setup. (b) Chip imaging device.




Figure S3. LAMP test for PDMS compatibility. (a) Left to right (bright field): Amplification result with mixing layer, amplification result with mixed layer-free and blank control. (b) Dark field image under UV excitation. (c) Electrophoresis of LAMP reaction. Lane M: DNA marker, lane 1: Blank control, lane 2: Amplification result with mixing layer. lane3: Amplification result mixed without mixed layer.




Figure S4. Specificity of LAMP reaction in the amplification tube. lane M: DNA marker, lane 1-3: *invA* primer test results. Lane 1: blank control, lane 2: *malB* template, lane 3: *invA* template. Lane 4-6: *malB* primer test results. Lane 4: blank control, lane 5: *invA* template, lane 6: *malB* template.