

Supporting Information

A complementary method with PFBBBr-derivatization based on GC-EI-MS platform for the simultaneous quantitation of short-, medium- and long -chain fatty acids in murine plasma and feces samples

Hong-Xu Zhou^{a,b#}, Qing Jiang^{a,b#}, Xin He^{a,b}, Xian Fu^{a,b}, Jun-Yan Liu^{a,b¹}

^aCNTTI of the Institute of Life Sciences & Anesthesia Department of the Second Affiliated Hospital, Chongqing Medical University, Chongqing, 400016, China

^bBasic Medicine Research and Innovation Center for Novel Target and Therapeutic Intervention, Ministry of Education, Chongqing, 400016, China

[#]These authors contributed equally to this work.

Table S1 The inter and intra-day of precision, accuracy for assays.

Table S2. The matrix effect, extraction recovery and stability for assays.

Table S3. The changes of FAs levels in different groups in feces and plasma

¹ To whom the correspondence should be addressed: Prof. Jun-Yan Liu, jyliu@cqmu.edu.cn, Tel; +86-23-6848 3587, Fax: +86-23-6366 4820

Table S1 The inter and intra-day of precision, accuracy for assays.

No.	Abbreviation	Inter-day assay (n=6)								Intra-day assay (n=18)							
		Precision (RSD%)				Accuracy (RE%)				Precision (RSD%)				Accuracy (RE%)			
		LLOQ	QCL	QCM	QCH	LLOQ	QCL	QCM	QCH	LLOQ	QCL	QCM	QCH	LLOQ	QCL	QCM	QCH
1	FA 2:0	1.42	0.78	2.61	0.89	106.55	99.67	95.70	90.73	4.39	6.79	3.23	4.16	109.86	103.33	101.73	94.04
2	FA 3:0	0.83	1.22	3.19	0.37	105.02	105.02	94.63	88.92	4.99	5.57	5.85	4.56	115.01	110.63	96.34	90.86
3	FA 4:1	1.03	1.16	3.47	0.50	108.61	99.05	94.28	91.69	5.39	5.47	7.51	5.42	117.03	111.16	93.76	92.28
4	FA 4:0	1.52	0.76	2.61	0.68	106.49	97.49	87.63	86.36	5.73	5.18	5.55	5.89	107.50	114.01	89.37	86.98
5	FA 5:1	2.53	3.80	3.37	0.56	115.53	105.53	94.41	93.43	4.91	11.52	7.02	3.91	116.78	114.22	94.65	96.49
6	FA 5:0	4.23	6.03	3.98	1.91	103.65	103.64	101.03	93.00	5.57	10.57	6.82	4.87	106.58	106.58	101.44	95.37
7	FA 6:0	2.37	1.51	9.08	1.21	112.50	102.51	86.07	91.60	17.07	5.69	13.23	5.70	117.37	110.13	88.32	92.44
8	FA 7:0	5.56	10.37	7.82	1.43	112.95	107.95	89.93	90.54	14.63	12.63	13.85	6.15	117.85	110.11	85.99	90.93
9	FA 8:0	5.99	7.71	7.29	0.93	103.69	113.69	107.66	111.72	11.88	13.54	11.53	4.24	111.59	101.59	114.75	112.74
10	FA 9:0	10.58	11.29	6.48	2.67	89.84	99.84	86.53	103.62	13.37	10.41	9.13	8.24	101.30	101.30	93.06	110.85
11	FA 10:0	6.74	0.91	1.61	1.76	100.34	109.74	110.13	113.11	7.95	8.91	11.69	4.96	114.14	102.14	113.45	113.36
12	FA 11:0	13.58	2.21	1.62	4.39	109.02	89.01	96.81	106.71	13.92	10.51	13.42	11.36	115.53	95.53	92.79	110.42
13	FA 12:0	11.32	1.33	11.96	2.56	109.47	99.47	87.49	104.52	14.81	14.07	14.73	11.81	112.35	92.35	85.51	106.89
14	FA 13:0	15.48	14.86	6.29	7.21	108.40	88.40	88.08	101.90	17.25	11.14	14.39	10.25	116.17	106.17	90.68	102.60
15	FA 14:1	14.58	2.48	0.79	8.72	109.16	90.17	95.37	92.44	16.99	5.25	7.58	18.85	110.37	114.38	104.09	90.59
16	FA 14:0	7.87	6.14	12.87	7.59	104.93	89.31	97.26	110.37	12.26	13.79	10.78	14.41	85.94	98.51	114.64	112.70
17	FA 15:1	9.90	3.68	8.76	8.06	99.33	93.31	90.57	85.18	17.82	14.25	10.64	14.86	118.69	98.71	87.03	89.05
18	FA 15:0	15.63	4.89	1.75	7.76	118.32	110.33	85.36	102.67	17.81	13.11	8.56	14.14	88.94	93.41	92.30	102.53
19	FA 16:1T	10.64	9.33	2.08	7.73	108.36	95.76	105.10	92.25	13.69	13.78	6.70	10.29	95.63	112.69	113.91	92.31
20	FA 16:1	5.03	10.79	3.07	7.80	116.35	109.47	97.67	95.91	19.55	10.92	8.49	10.36	117.63	112.59	108.11	95.52
21	FA 16:0	4.38	12.75	9.33	4.02	83.36	89.23	98.89	106.18	9.16	7.60	11.29	10.54	82.93	86.46	88.60	111.34
22	FA 17:1T	10.66	5.01	0.42	5.59	95.36	89.15	113.23	95.71	12.56	12.60	5.12	13.08	99.88	97.62	108.88	96.83

23	FA 17:0	11.04	3.41	4.14	7.19	102.33	88.12	112.31	98.85	13.25	7.69	6.24	10.41	102.44	110.88	110.85	98.82
24	FA 18:3N6	12.90	8.02	2.47	3.78	103.65	99.59	104.29	89.96	16.30	14.61	14.03	6.31	118.43	111.99	114.44	87.68
25	FA 18:2N6	13.84	13.08	1.74	4.14	112.38	102.51	86.39	89.61	15.61	12.79	6.16	8.53	112.56	108.95	85.26	86.85
26	FA 18:1N9C	4.30	3.51	1.92	5.90	89.36	88.25	87.03	89.62	9.77	11.63	5.44	12.84	110.89	96.32	85.23	85.16
27	FA 18:3N3	12.62	9.38	2.35	5.60	116.29	114.22	86.48	86.66	15.33	11.17	10.78	11.01	109.68	114.32	86.42	90.98
28	FA 18:1N9T	10.18	9.16	2.88	6.48	109.94	105.69	99.49	87.06	18.15	12.07	11.26	12.23	99.59	96.59	94.07	93.80
29	FA 18:1N7	11.88	9.43	2.86	5.48	118.49	110.06	114.95	89.38	14.51	13.26	8.17	12.67	107.59	95.33	113.43	85.28
30	FA 18:0	8.63	13.04	11.27	6.43	109.58	88.69	112.48	105.85	9.88	10.40	12.44	13.84	110.39	96.32	112.96	109.79
31	FA 19:1N9T	13.64	6.67	2.19	3.91	106.89	112.63	102.27	88.25	17.71	9.57	3.73	8.37	110.95	112.17	102.89	89.61
32	FA 19:0	10.45	6.11	1.32	4.69	116.69	95.03	103.56	86.65	14.52	7.32	4.80	10.66	89.65	110.87	102.30	86.06
33	FA 20:4N6	11.53	10.29	4.35	3.04	112.45	96.25	102.39	86.71	13.93	14.43	5.61	4.05	109.75	103.86	102.70	87.09
34	FA 20:5N3	9.87	10.98	5.91	5.25	104.75	92.51	89.17	89.14	14.25	10.47	7.78	3.75	100.85	90.25	86.15	91.18
35	FA 20:1T	10.07	12.31	9.44	3.46	112.59	105.57	111.59	93.13	12.91	14.47	10.44	7.06	109.59	114.18	114.32	91.88
36	FA 20:1	8.09	12.71	7.67	4.05	116.35	90.23	114.59	88.19	13.95	14.39	10.52	8.38	115.84	95.11	112.93	85.79
37	FA 20:2	10.34	13.35	4.27	5.76	98.56	105.05	107.85	87.66	11.50	14.11	9.56	9.49	102.59	114.94	108.68	85.06
38	FA 20:0	9.01	1.79	3.49	5.40	94.55	101.84	95.91	91.54	10.16	7.10	6.16	12.26	112.86	110.57	93.73	87.09
39	FA 21:0	5.89	5.00	2.93	3.44	108.55	98.45	109.58	86.61	8.74	7.59	13.75	8.39	109.63	93.46	106.03	87.70
40	FA 22:0	14.61	12.69	9.92	13.32	81.29	95.82	110.14	97.31	19.30	12.00	13.99	13.01	107.59	94.25	100.18	98.02
41	FA 22:6N3	12.29	11.21	3.96	1.82	80.99	112.81	95.23	87.12	16.89	11.06	7.07	4.74	100.29	94.22	104.76	86.42
42	FA 22:2	10.91	5.89	0.82	3.33	109.28	104.32	98.99	87.96	11.80	6.50	3.97	8.42	114.87	112.63	98.87	86.18
43	FA 23:0	9.05	6.71	2.24	3.25	99.57	104.09	101.48	88.23	9.96	9.08	4.74	8.23	118.52	110.58	100.68	88.71
44	FA 24:0	10.52	6.05	2.16	2.32	96.59	101.89	93.03	86.74	12.73	5.41	3.18	7.39	110.35	100.15	94.06	89.76

Table S2. The matrix effect, extraction recovery and stability for assays.

NO.	Abbreviation	Extraction recovery (%)			Matrix effect (%)			Stability (%)	
		QCL	QCM	QCH	QCL	QCM	QCH	24 h	48 h
1	FA 2:0	89.53	118.79	117.73	89.63	86.46	87.25	92.48	92.64
2	FA 3:0	86.14	117.49	88.00	86.97	85.44	89.22	91.68	92.11
3	FA 4:1	83.19	94.82	109.12	85.63	91.80	99.01	90.26	90.80
4	FA 4:0	89.01	87.10	103.58	86.51	88.09	94.52	89.34	89.74
5	FA 5:1	100.41	101.79	108.25	102.78	97.81	98.21	93.72	94.26
6	FA 5:0	85.56	87.41	108.01	88.12	86.82	97.91	93.06	93.82
7	FA 6:0	83.42	104.42	111.56	90.02	101.93	101.36	92.08	93.16
8	FA 7:0	89.03	103.93	112.37	89.30	100.44	101.99	91.88	92.98
9	FA 8:0	81.41	104.96	108.42	85.98	88.81	86.55	93.71	94.40
10	FA 9:0	85.08	104.16	116.12	88.98	88.01	92.77	89.45	91.32
11	FA 10:0	94.66	115.81	109.47	110.48	99.08	87.57	93.13	93.85
12	FA 11:0	98.98	115.29	120.71	100.34	104.59	98.06	86.67	89.34
13	FA 12:0	88.21	106.13	114.81	86.32	89.86	99.62	84.65	87.56
14	FA 13:0	88.27	120.39	119.25	91.17	101.44	104.50	82.41	85.88
15	FA 14:1	92.61	108.29	117.36	95.53	107.60	114.11	77.26	81.57
16	FA 14:0	96.20	86.53	118.02	87.09	90.35	105.12	84.96	87.82
17	FA 15:1	105.21	110.58	118.69	113.66	114.89	110.90	73.11	77.00
18	FA 15:0	100.15	103.06	113.95	105.50	114.34	113.44	80.19	84.00
19	FA 16:1T	79.09	116.84	117.01	102.29	112.35	103.61	80.51	83.08
20	FA 16:1	103.45	107.50	116.55	92.44	97.31	110.60	79.24	83.45
21	FA 16:0	93.28	102.97	114.98	117.14	111.73	108.38	86.48	88.84
22	FA 17:1T	83.83	106.38	97.85	86.24	89.83	113.42	81.30	84.03
23	FA 17:0	100.57	120.23	153.62	112.23	113.44	104.83	79.90	83.38

24	FA 18:3N6	96.43	97.13	118.13	91.07	89.96	112.23	90.67	93.49
25	FA 18:2N6	110.72	117.27	118.16	105.48	103.11	108.36	87.34	88.65
26	FA 18:1N9C	88.32	112.17	119.53	108.73	107.12	112.91	81.64	84.01
27	FA 18:3N3	101.35	109.29	112.15	104.82	106.57	100.45	87.37	89.37
28	FA 18:1N9T	110.92	115.26	118.64	98.14	86.42	100.14	83.07	85.79
29	FA 18:1N7	108.90	104.85	117.55	103.36	86.25	112.25	83.12	85.93
30	FA 18:0	97.58	83.85	116.84	99.37	90.30	107.22	80.42	83.95
31	FA 19:1N9T	108.78	77.35	117.46	103.17	85.19	112.88	87.63	89.27
32	FA 19:0	106.40	110.91	119.58	110.35	110.23	102.61	84.33	86.20
33	FA 20:4N6	104.03	107.26	113.31	102.77	110.96	114.34	92.77	93.35
34	FA 20:5N3	87.05	107.80	112.61	119.55	113.28	112.89	98.00	89.99
35	FA 20:1T	105.03	96.00	113.45	104.58	106.40	113.18	92.01	93.53
36	FA 20:1	81.35	110.11	117.86	113.39	104.11	109.21	88.46	89.59
37	FA 20:2	91.96	88.95	116.80	89.43	108.29	111.60	88.92	87.94
38	FA 20:0	113.30	103.61	114.95	110.81	112.52	113.22	83.67	85.81
39	FA 21:0	88.24	90.58	113.13	89.37	82.58	102.18	86.94	88.08
40	FA 22:0	107.22	117.47	119.23	88.31	113.57	111.78	95.13	95.17
41	FA 22:6N3	80.32	91.08	108.73	103.25	109.85	114.21	93.45	93.39
42	FA 22:2	109.18	108.45	115.82	106.52	111.06	107.09	87.04	88.36
43	FA 23:0	104.97	114.50	112.60	102.09	106.69	103.67	87.03	88.34
44	FA 24:0	113.36	114.41	111.67	113.28	107.39	102.63	88.70	89.11

Table S3. The changes of FAs levels in different groups in feces and plasma

No.	Abbreviation	Control group (mean \pm SD)	DON group (mean \pm SD)	Control group versus DON group		
				VIP value ^a	<i>p</i> -value ^b	Fold change
1	FA 2:0_F	83.21 \pm 46.67	141.30 \pm 51.86	0.8353	0.0260	0.59
2	FA 4:1_F	1.80 \pm 1.19	0.95 \pm 0.48	0.8729	0.0485	1.89
3	FA 4:0_F	6.79 \pm 5.84	17.26 \pm 9.41	0.8011	0.0062	0.39
4	FA 10:0_F	3.53 \pm 2.07	11.14 \pm 4.18	1.1898	0.0005	0.32
5	FA 14:0_F	2.75 \pm 0.62	8.66 \pm 2.24	1.5209	<0.0001	0.32
6	FA 15:0_F	3.96 \pm 1.63	10.93 \pm 3.61	0.4692	<0.0001	0.36
7	FA 16:1T_F	1.72 \pm 0.41	3.12 \pm 0.81	1.3185	0.0004	0.55
8	FA 16:0_F	71.68 \pm 20.05	171.60 \pm 32.69	1.4937	<0.0001	0.42
9	FA 17:0_F	3.36 \pm 0.72	7.63 \pm 0.73	1.5963	<0.0001	0.44
10	FA 18:1N9C_F	27.78 \pm 10.93	57.83 \pm 16.18	1.3110	0.0003	0.48
11	FA 18:1N9T_F	4.23 \pm 1.37	11.26 \pm 4.62	1.3051	0.0004	0.38
12	FA 18:0_F	73.10 \pm 12.98	129.43 \pm 15.09	1.5023	<0.0001	0.56
13	FA 19:0_F	1.65 \pm 0.38	2.31 \pm 0.32	1.1276	0.0010	0.71
14	FA 20:0_F	18.00 \pm 4.51	21.64 \pm 2.75	0.7521	0.0569	0.83
15	FA 23:0_F	6.07 \pm 1.64	8.76 \pm 0.87	1.1864	0.0002	0.69
16	FA 10:0_P	36.24 \pm 16.61	51.17 \pm 11.21	1.4051	0.0211	0.71
17	FA 12:0_P	11.50 \pm 10.39	18.64 \pm 9.50	1.0834	0.0894	0.62
18	FA 14:1_P	0.89 \pm 0.23	0.69 \pm 0.19	1.2145	0.0836	1.29
19	FA 16:1_P	1.81 \pm 2.08	0.25 \pm 0.36	1.4599	0.0301	7.34
20	FA 17:0_P	0.16 \pm 0.05	0.23 \pm 0.10	1.3034	0.0340	0.67
21	FA 20:4N6_P	2.89 \pm 0.75	2.02 \pm 0.55	1.5857	0.0095	1.43

a The VIP value was from the OPLS mode.

b The *p*-value was calculated from the Student's *t*-test.