

## Supporting Information

# Sustainable pathway to furanics from biomass via heterogeneous organo-catalysis

Sanny Verma<sup>a†</sup>, R. B. Nasir Baig<sup>a†</sup>, Mallikarjuna N. Nadagouda,<sup>b</sup> Christophe Len<sup>c</sup> and Rajender

S. Varma<sup>a\*</sup>

<sup>a</sup>Sustainable Technology Division, National Risk Management Research Laboratory, U. S. Environmental Protection Agency, MS 443, Cincinnati, Ohio 45268, USA. Fax: 513- 569-7677; Tel: 513-487-2701. E-mail: [varma.rajender@epa.gov](mailto:varma.rajender@epa.gov)

<sup>b</sup>WQMB, WSWRD, National Risk Management Research Laboratory, U. S. Environmental Protection Agency, Cincinnati, Ohio 45268, USA

<sup>c</sup>Sorbonne Universités, Université de Technologie de Compiègne, Compiègne, France

1. Synthesis of g-CN and Sg-CN
  - (a) Synthesis of g-CN
  - (b) Synthesis of Sg-CN
2. General procedure for the synthesis of carbohydrates to value added chemicals
3. General procedure for the synthesis of levulinic acid from glucose
4. Recycling of Sg-CN catalyst (**S1**)
5. Higher magnification SEM image of Sg-CN catalyst (**S2**)
6. TEM image (different view) of Sg-CN catalyst (**S3**)
7. XRD of Sg-CN catalyst (**S4**)
8. N<sub>2</sub> sorption isotherms of Sg-CN (**S5**)
9. Distribution of pore diameter of Sg-CN (**S6**)
10. N<sub>2</sub> sorption isotherms of g-CN (**S7**)
11. Distribution of pore diameter of g-CN (**S8**)
12. EDX of Sg-CN catalyst (**S9**)
13. GC-MS data of the product
14. <sup>1</sup>H and <sup>13</sup>C NMR of the product

## **1. Synthesis of g-CN and Sg-CN catalyst**

### **a) Synthesis of g-CN**

The pure urea obtained from Aldrich was calcinated at 500 °C for 2 hours in a closed furnace. A pale-yellow solid of pure graphitic carbon nitride (g-CN) was obtained and used without any purification.

### **b) Synthesis of Sg-CN catalyst**

Graphitic carbon nitride, g-CN (1.0 g) and dichloromethane (50 mL) were taken in a round bottom flask. Chlorosulfonic acid (0.5 mL) was added to the reaction mixture over the period of 10 min under continuous stirring and then stirring was continued for 3 hours. The ensuing white solid was filtered off, washed with water, methanol and dried under vacuum at 50 °C. The Sg-CN catalyst was characterized by transmission electron microscopy (TEM), scanning electron microscopy (SEM), X-ray diffraction (XRD), Fourier transform infrared spectroscopy (FTIR), Energy-dispersive X-ray spectroscopy (EDX), X-ray photoelectron spectroscopy (XPS) and thermogravimetric analysis (TGA).

## **2. General procedure for the synthesis of carbohydrates to value added chemicals**


A reaction tube equipped with a stir bar was charged with the desired amount of carbohydrate (2 mmol), catalyst (25 mg of Sg-CN and 5.0 mg of KBr, if required) and water (2 mL). The glass tube was then heated in oil bath at 100 °C. After the completion of the reaction, the reaction temperature was brought down to room temperature, the catalyst was recovered *via* filtration/centrifugation and the product was isolated using solvent extraction using ethyl acetate, dried over sodium sulfate, concentrated, and characterized.

### **3. General procedure for the synthesis of levulinic acid from glucose**

A reaction tube equipped with a stir bar was charged with the desired amount of glucose (2 mmol), catalyst (25 mg), and water (2 mL). The glass tube was sealed and placed in an oil bath and heated at 150 °C over a period of 8 hours. After completion of the reaction, the reaction temperature was allowed to come at room temperature. The catalyst was recovered using a centrifuge and the product was isolated using solvent extraction.

#### 4. Recycling of the catalyst

After the completion of each reaction, the Sg-CN catalyst was recovered using a centrifuge,




**S1.** Recycling of the catalyst

washed with water followed by acetone, and dried under vacuum and used for the fresh set


of reactants. It was observed that the catalyst remains active even 3<sup>rd</sup> cycle of the reaction.

5. Higher magnification SEM image of Sg-CN catalyst




S2. Higher magnification SEM image of Sg-CN catalyst

6. TEM image (different view) of Sg-CN catalyst.




S3. TEM image (different view) of Sg-CN catalyst.

## 7. XRD of Sg-CN catalyst



S4. XRD of Sg-CN catalyst


## 8. N<sub>2</sub> sorption isotherms of Sg-CN



S5. N<sub>2</sub> sorption isotherms of Sg-CN




## 9. Distribution of pore diameter of Sg-CN



S6. Distribution of pore diameter of Sg-CN


## 10. N<sub>2</sub> sorption isotherms of g-CN



S7. N<sub>2</sub> sorption isotherms of g-CN




### 11. Distribution of pore diameter of g-CN



S8. Distribution of pore diameter of g-CN



## 12. EDX of Sg-CN catalyst





S9. EDX of Sg-CN catalyst

### 13. GC-MASS data of the product



File :C:\Sanny\Data\NS 125 GH.D  
Operator : kt  
Acquired : 20 Oct 2015 12:43 pm using AcqMethod SAMPLES- 20 MIN\_20151015\_A.M  
Instrument : Instrument #1  
Sample Name: NS 125 GH  
Misc Info :  
Vial Number: 12



File :C:\Sanny\Data\NS 128.D  
Operator : kt  
Acquired : 26 Oct 2015 3:53 pm using AcqMethod SAMPLES- 20 MIN\_20151015\_A.M  
Instrument : Instrument #1  
Sample Name: NS 128  
Misc Info :  
Vial Number: 15





File :C:\Sanny\Data\ns 128 gyt.D  
Operator : kt  
Acquired : 3 Nov 2015 3:31 pm using AcqMethod SAMPLES- 20 MIN\_20151015\_A.M  
Instrument : Instrument #1  
Sample Name: ns 128 gyt  
Misc Info :  
Vial Number: 6








File :C:\Sanny\Data\NS 135.D  
Operator : kt  
Acquired : 25 Nov 2015 4:13 pm using AcqMethod SAMPLES- 20 MIN\_20151015\_A.M  
Instrument : Instrument #1  
Sample Name: NS 135  
Misc Info :  
Vial Number: 13




File :C:\Sanny\Data\NS 164 (35).D  
Operator : kt  
Acquired : 25 Jan 2016 11:20 am using AcqMethod SAMPLES- 20 MIN\_20151015\_A.M  
Instrument : Instrument #1  
Sample Name: NS 164 (35)  
Misc Info :  
Vial Number: 100




## 14. $^1\text{H}$ and $^{13}\text{C}$ NMR of the product




NS 124 13CNMR



NS 125








NS 128 F

— 9.885


— 7.349




NS 128 13CNMR










NS 135\_1HNMR

7.870




2.779  
2.776  
2.757  
2.734  
2.647  
2.644  
2.627  
2.623  
2.621  
2.604  
2.600




```
NAME      NS 135 1HNMR
EXPNO     1
PROCNO    1
Date_     2015102
Time      15.48
INSTRUM   spect
PROBHD    5 mm BBO BB-1H
PULPROG   zgpg30
TD         32768
SOLVENT   CDCl3
NS         16
DS         2
SWH        6188.119 Hz
FIDRES    0.198846 Hz
AQ         2.6477852 sec
RG         181
DM         80.800 usec
DE         6.50 usec
TE         303.0 K
D1         1.00000000 sec
TDO        1
```

```
===== CHANNEL f1 =====
NUC1       1H
P1         13.00 usec
PL1        3.20 dB
PL1W       12.02264404 W
SFO1       300.1318534 MHz
SI         32768
SF         300.1300000 MHz
WDW        EM
SSB        0
LB         0.30 Hz
GB         0
PC         1.00
```



NS 135\_13CNMR

207.18  
177.92



37.69  
29.67  
27.75



NAME NS 135\_13CNMR  
EXPNO 2  
PROCNO 1  
Date\_ 20151102  
Time 21.09  
INSTRUM spect  
PROBHD 5 mm BBO BB-1H  
PULPROG zgpg  
TD 32768  
SOLVENT CDCl3  
NS 5000  
DS 4  
SWH 17985.611 Hz  
FIDRES 0.548877 Hz  
AQ 0.9110282 sec  
RG 8192  
DW 27.800 usec  
DE 6.50 usec  
TE 303.0 K  
D1 2.0000000 sec  
D11 0.0300000 sec  
TD0 1

----- CHANNEL f1 -----  
NUC1 13C  
P1 10.00 usec  
PL1 1.80 dB  
PLLW 49.78760910 W  
SFO1 75.4752953 MHz

----- CHANNEL f2 -----  
CPDPRG2 waltz16  
NUC2 1H  
PCPD2 80.00 usec  
PL2 3.20 dB  
PLL2 18.98 dB  
PL2W 12.02264604 W  
PLL2W 0.31768745 W  
SFO2 300.1312005 MHz  
S1 32768  
SFO 75.4677490 MHz  
WDW EM  
SSB C  
TE 1.00 Hz  
CB C  
PC 1.40

