

Supporting Information

Analysis of trace elements in uranium by inductively coupled plasma optical emission spectroscopy, design of experiments, and partial least squares regression

Luke R. Sadergaski, Benjamin T. Manard, and Hunter B. Andrews,*

Oak Ridge National Laboratory, 1 Bethel Valley Rd., Oak Ridge, TN, 37980, USA

*Corresponding author: sadergaskilr@ornl.gov

This Supporting Information contains:

Pages: 4

Figures: 2

Tables: 4

Table S1. Elements in one trace 100 $\mu\text{g mL}^{-1}$ standard.

Al	Cr	Fe	K	Tl
As	Co	La	Pr	Th
Ba	Cu	Pb	Re	Tm
Be	Dy	Li	Rb	U
Bi	Er	Lu	Sm	V
B	Eu	Mg	Sc	Yb
Cd	Gd	Mn	Se	Y
Ca	Ga	Nd	Na	Zn
Ce	Ho	Ni	Sr	
Cs	In	P	Tb	

Table S2. Elements in the second trace 100 $\mu\text{g mL}^{-1}$ standard.

Sb	Ti
Ge	W
Hf	Zr
Mo	
Nb	
Si	
Ag	
Ta	
Te	
Sn	

Table S3. Additional samples included in the validation set (*i.e.*, quality controls).

Sample	Uranium ($\mu\text{g mL}^{-1}$)	Trace ($\mu\text{g mL}^{-1}$)	($\mu\text{g/g}$)
1	1000	0.02	20
2	1000	0.05	50
3	1000	0.1	100
4	1000	0.25	250
5	1000	0.5	500
6	1000	1.5	1500




Figure S1. Optical emission spectra of trace elements from 0.02 to 1.5 $\mu\text{g mL}^{-1}$ in 1000 $\mu\text{g mL}^{-1}$ U for (a) Mn 259.37 nm, (b) Mn 257.61 nm, (c) Fe 238.20 nm, and (d) Fe 259.94 nm. Compared to a standard (1 $\mu\text{g mL}^{-1}$ multielement) and a 1000 $\mu\text{g mL}^{-1}$ U sample.




Figure S2. Optical emission spectra of Nb 309.42 nm 0.02 to 1.5 $\mu\text{g mL}^{-1}$ in 1000 $\mu\text{g mL}^{-1}$ U. Compared to a standard (1 $\mu\text{g mL}^{-1}$ multielement) and a 1000 $\mu\text{g mL}^{-1}$ U sample.

Table S4. Quality control samples to determine RMSE.

Sample	Trace ($\mu\text{g mL}^{-1}$)	U ($\mu\text{g mL}^{-1}$)	$\mu\text{g trace/g U}$
1	1.164	528.3	2208
2	1.454	762.4	1911
3	0.543	107.9	5055
4	0.776	344.5	2258
5	0.020	963.9	20
6	0.052	962.8	54
7	0.099	966.0	102
8	0.241	965.3	250
9	0.482	967.3	498
10	1.569	975.6	1608