

Supporting Information

Development of biomass waste-based carbon quantum dots and their potential application as non-toxic bioimaging agents

Norhidayah Abu^{1,2*}, Shanmugavel Chinnathambi^{3*}, Mahima Kumar³, Fatemeh Etezadi³,
Noremylia Mohd Bakhori², Zuhana Ahmad Zubir², Shahrul Nizam Md Salleh², Rafidah Hanim
Shueb¹, Subramani Karthikeyan⁴, Vaijayanthi Thangavel³, Jaafar Abdullah^{5*} and Ganesh N.
Pandian^{3*}.

¹Department of Medical Microbiology & Parasitology, School of Medical Sciences, Universiti Sains Malaysia, Health Campus, Kubang Kerian 16150, Kelantan, Malaysia


²Advanced Materials Research Centre (AMREC), SIRIM Berhad, Lot 34, Jalan Hi-Tech 2/3, Kulim, Hi-Tech Park, 09000 Kulim, Malaysia

³Institute for Integrated Cell-Material Sciences, Institute for Advanced Study, Kyoto University, Kyoto, Japan

⁴Centre for Healthcare Advancement, Innovation and Research, Vellore Institute of Technology, Chennai, 600 127, India


⁵Department of Chemistry, Faculty of Science, Universiti Putra Malaysia, UPM Serdang, Serdang 43400, Selangor, Malaysia

C.A. E-mail: nhidayah@sirim.my; chinnathambi.shanmugavel.8s@kyoto-u.ac.jp;
jaafar@upm.edu.my; namasivayam.ganeshpandian.5z@kyoto-u.ac.jp




Sample (n=3)	Quantum Yield (%)	Abs	Peak Wavelength	Peak Count	Peak FWHM
Reference (Series 1)	-----	-----	371.88	65.27	141.66
Palm shell CDs 1 (Series 2)	2.4	0.914	385.5	81.26	74.74
Palm shell CDs 2 (Series 3)	2.4	0.914	390.04	83.2	69.04
Palm shell CDs 3 (Series 4)	2.3	0.914	388.53	82.88	71.34

Fig.S1 Palm shell carbon quantum dots and their absolute quantum yield measurements: Figure and its related parameters.



Sample (n=3)	Quantum Yield (%)	Abs	Peak Wavelength	Peak Count	Peak FWHM
Reference (Series 1)	---	---	392.31	36.51	423.69
Oyster shell CDs 1 (Series 2)	1.4	0.967	393.07	41.72	89.03
Oyster shell CDs 2 (Series 3)	1.3	0.967	403.65	40.09	90.84
Oyster shell CDs 3 (Series 4)	1.4	0.967	392.31	41.03	88.79

Fig. S2 Oyster shell carbon quantum dots and their absolute quantum yield measurements: Figure and its related parameters.



Sample (n=3)	Quantum Yield (%)	Abs	Peak Wavelength	Peak Count	Peak FWHM
Reference (Series 1)	---	---	406.67	26.74	165.51
CA CDs 1 (Series 2)	22.1	0.986	461.72	381.47	102.29
CA CDs 2 (Series 3)	22.1	0.986	466.24	381.01	102.33
CA CDs 3 (Series 4)	22.1	0.986	465.49	381.46	103.04

Fig. S3 Citric acid carbon quantum dots and their absolute quantum yield measurements: Figure and its related parameters.




Fig. S4 Effect of freeze-dry process vs. absorbance of carbon dots




Fig. S5 The excitation spectrum of CPKS carbon quantum dots (yellow line) and its emission spectra with various wavelengths (excitation at 300 nm, 315nm, 325nm)




Fig.S6 The excitation spectrum of oyster carbon quantum dots (yellow line) and its emission spectra with various wavelengths (excitation at 300 nm, 325nm, 350nm)




Fig. S7 The excitation spectrum of citric acid carbon quantum dots (yellow line) and its emission spectra with various wavelengths (excitation at 325 nm, 350nm, 375nm)




Fig. S8 XPS full-survey spectrum of CQDs from CPKS (a); high-resolution XPS of the C1s (b), N1s (c), and O1s (d) spectra.




Fig. S9 XPS full-survey spectrum of CQDs from Oyster shell (a); high-resolution XPS of the C1s (b), N1s (c), and O1s (d) spectra.




Fig. S10 XPS full-survey spectrum of CQDs from citric acid (a); high-resolution XPS of the C1s (b), N1s (c), and O1s (d) spectra.

Table.S1 C1s, N1s, and O1s and the relative atomic percentages of these peaks

Type of CQDs	Binding energy (eV)							Element content (%)		
	C-C, C-H	C-N	C-N, C-O	C=O, OH-C-H	O=C	O-C	R-NH-R	C	O	N
CPKS	284.68	285.87	287.90	289.61	530.95	532.35	399.3	58.39	32.35	9.26
Oyster shell	284.63	285.78	288.08		531.34	533.29	399.11	65.54	30.94	3.52
Citric acid	284.79	286.73			530.84	532.03	399.92	76.87	12.35	10.78




Fig. S11 Zeta potential measurement of synthesized QDs from three different sources.




Fig. S12 Mitochondrial localization using CPKS, Oyster, and Citric acid-based CQDs with HeLa cells.




Fig. S13 Medaka egg (day-9) with increasing concentrations of CPKS - QDs 250, 500, 750 and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S14 Medaka egg (day-11) with increasing concentrations of CPKS - QDs 250, 500, 750 and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S15 Medaka egg (day-14) with increasing concentrations of CPKS - QDs 250, 500, 750 and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S16 Medaka egg (day-9) with increasing concentrations of **Oyster - QDs** 250, 500, 750 and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S17 Medaka egg (day-11) with increasing concentrations of **Oyster - QDs** 250, 500, 750, and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S18 Medaka egg (day-14) with increasing concentrations of **Oyster - QDs** 250, 500, 750 and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S19 Medaka egg (day-9) with increasing concentrations of CA-QDs 250, 500, 750 and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S20 Medaka egg (day-11) with increasing concentrations of CA-QDs 250, 500, 750, and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).




Fig. S21 Medaka egg (day-14) with increasing concentrations of CA-QDs 250, 500, 750 and 1000 $\mu\text{g mL}^{-1}$. Three eggs were used for each micelle concentration \pm S.E. ($n = 3$).