

## Supplementary Material


### Degradation and detection of organophosphorus pesticides based on peptides and MXene-peptide composite materials

Qiuying Wang<sup>1a</sup>, Ruiqing Jiu<sup>1a</sup>, Yunyao Wang<sup>a</sup>, Zongda Li<sup>a</sup>, Jianan Chen<sup>a</sup>, Haochi Liu<sup>a</sup>, Jifeng Liu<sup>\*a</sup>, Jia Cao<sup>b</sup>


*a: State Key Laboratory of Food Nutrition and Safety, Key Laboratory of Food Quality and Healthy of Tianjin, College of Food Science and Engineering, Tianjin University of Science and Technology, Tianjin 300457, PR China*

*b: Tianjin Vocational College of Bioengineering*


*\*Corresponding author: Jifeng Liu, liujifeng111@gmail.com*




**Scheme S1** Mechanism of OPs hydrolysis catalyzed by enzyme. X<sub>1</sub>, X<sub>2</sub>, and X<sub>3</sub> represent acids, bases, and nucleophiles, respectively.



**Fig. S1** Mass spectrum of peptides.



**Fig. S2** UV-vis absorbance spectra of PNP ([PNP] = 0.1 mM).



**Fig. S3** Geometric configuration of peptides (carbon atoms: gray; nitrogen atoms: blue; oxygen atoms: red; hydrogen atoms: white).




Fig. S4. The ATR-FTIR spectra of MXene, MXene-CA and MXene-peptide.




Fig. S5 The oxidation of PNP.




Fig. S6 The standard curve of gas chromatography-mass spectrometry of OPs.

**Table S1** Compared to other detection methods

Determinative method	LOD	Linear range	RSD	Detection time
MXene-peptide electrochemical detection (this work)	0.15 $\mu\text{M}$	1-100 $\mu\text{M}$	4.71%	15 min
Roll-to-Roll manufactured sensors <sup>1</sup>	1 $\mu\text{M}$	1-20 $\mu\text{M}$	<10%	10 min
Fluorescent peptide probes <sup>2</sup>	0.6 $\mu\text{M}$	1-100 $\mu\text{M}$	<2%	15 min
Colourimetric sensing platform <sup>3</sup>	0.39 $\mu\text{M}$	0.5-50 $\mu\text{M}$	-	2 min
Electrochemical detection based on amino acids <sup>4</sup>	0.24 $\mu\text{M}$	0.5-100 $\mu\text{M}$	2.3%	15 min
Multi-enzyme/CNT biosensor <sup>5</sup>	0.50 $\mu\text{M}$	0.5-40 $\mu\text{M}$	3.6%	>15 min

**Table S2** Analysis of OPs in real samples.

Sample	OPs	Additive amount ( $\mu\text{M}$ )	Initial value	GC-MS		Electrochemical analysis			
				Average value ( $\mu\text{M}$ )	Recovery(%)	RS D (%) n=3	Average value ( $\mu\text{M}$ )	Recovery(%)	RS D (%) n=3
Oilseed rape	methyl paraoxon	5.00	0	4.47	85.59	3.66	4.10	82.00	4.71
		20.00	0	17.68	87.58	3.45	17.66	88.30	4.43
		40.00	0	32.34	77.72	4.09	35.64	89.10	4.30
Bokchoy	methyl paraoxon	5.00	0	4.08	74.65	4.32	3.84	76.80	4.15
		20.00	0	17.28	85.41	3.93	16.87	84.35	3.54
		40.00	0	32.51	78.17	4.21	35.38	88.45	2.95

## References

1. A. M. Ulloa, N. Glassmaker, M. R. Oduncu, P. Xu, A. Wei, M. Cakmak and L. Stanciu, *ACS Applied Materials & Interfaces*, 2021, **13**, 35961-35971.
2. J. Wang, J. Zhang, J. Wang, G. Fang, J. Liu and S. Wang, *Journal of Hazardous Materials*, 2020, **389**, 122074.
3. C.-L. Yang, L.-H. Yu, Y.-H. Pang and X.-F. Shen, *Analytica Chimica Acta*, 2024, **1286**, 342045.
4. Y. Liu, X. Cao, Z. Liu, L. Sun, G. Fang, J. Liu and S. Wang, *Analyst*, 2021, **145**, 8068-8076.
5. Y. Zhang, M. A. Arugula, M. Wales, J. Wild and A. L. Simonian, *Biosensors and Bioelectronics*, 2015, **67**, 287-295.