

**Supplementary Table S3 UPLC elution program applied to the separation
of *O*-glycans**

Time (min)	Flow rate (mL/min),	Concentration of ammonium formate (%)	Concentration of acetonitrile (%)
0.00	0.50	5.00	95.00
6.00	0.50	22.00	78.00
44.50	0.50	44.10	55.90
46.50	0.25	100.00	0.00
48.50	0.25	100.00	0.00
52.00	0.25	5.00	95.00
55.00	0.50	5.00	95.00
57.00	0.50	5.00	95.00