

Supporting information for

Green light-responsive photoelectrochemical sensing nanoplatform based on copper cobaltite nanorods for ultrasensitive detection of furazolidone antibiotic residue in food samples

Nguyen Ngoc Huyen^{1,*}, Ngo Xuan Dinh¹, Vu Thi Huong Mai¹, Phung Thi Lan Huong¹,

Pham Duc Thang², Anh-Tuan Le^{**1,3}

(1) Phenikaa University Nano Institute (PHENA), PHENIKAA University, Hanoi 12116, Vietnam.

(2) Faculty of Physics, VNU Hanoi University of Science, Vietnam National University, Hanoi 11416, Vietnam

(3) Faculty of Materials Science and Engineering, PHENIKAA University, Hanoi 12116, Vietnam.

Corresponding authors:

*huyen.nguyennngoc@phenikaa-uni.edu.vn (N. N. Huyen)

**tuan.leanh@phenikaa-uni.edu.vn (A.-T. Le)

Content

S1. Supplemental Results and Discussions




Figure S1. (a - b) Experimental and fitted Nyquist plots of impedance spectra corresponding in the frequency range from 0.01 Hz to 50 kHz of bare SPE, and CCO/SPE. (c -d) The Randles equivalent circuit used for fitting the data.




Figure S2. (a) DPV voltammograms of CCO/SPE in 200 μ M FZD at different pH values from 3.0 to 11.0; (b) the corresponding calibration plots of electro-reduction peak current and E_{pc} of FZD against pH values with error bars, (c) DPV curves of CCO/SPE in 200 μ M FZD at various accumulation time values, and (d) current response with accumulation time values from DPV measurements.

Table S1. A comparison of the electrochemical performance of different modified electrodes toward FZD detection.

Electrode	Technique	Linear range (μM)	LOD (μM)	Sensitivity ($\mu\text{A} \mu\text{M}^{-1} \text{cm}^{-2}$)	Real samples	Ref.
ZnFe ₂ O ₄ /SPE	DPV	0.5 - 75	0.61	0.56	Milk	¹
ZnO/ZnFe ₂ O ₄ /SPE	DPV	0.5 - 75	0.41	0.70		
ZnO/ZnFe ₂ O ₄ /SPE	DPV	1 - 100	0.65	0.76	Pork and shrimp	²
GO/SPE	LSV	1 - 100	0.19	1.03	-	³
Ag@Fe ₃ O ₄ /SPE	DPV	0.5 - 100	0.24	1.36	Honey and milk	⁴
NiFe ₂ O ₄ /rGO/GCE	DPV	0.1 - 150	0.05	-	Plasma, urine	⁵
CuCo ₂ O ₄ /SPE (Without green light illumination)	DPV	0.5 - 200	0.09	0.93	-	This work
CuCo ₂ O ₄ /SPE (With green light illumination)	DPV	0.25 -200	0.03	1.11	Honey and milk	

LOD: Limit of detection; **SPE:** Screen-printed electrode; **DPV:** Differential pulse voltammetry; **GO:** Graphene oxide; **LSV:** Linear sweep voltammetry; **rGO:** Reduced graphene oxide; **GCE:** Glassy carbon electrode;

Reference

1. N. N. Huyen, N. T. Anh, T. L. H. Phung, N. X. Dinh, N. T. Vinh, T. T. Loan, D. L. Vu and A.-T. Le, Boosting the Selective Electrochemical Signals for Simultaneous Determination of Chloramphenicol and Furazolidone in Food Samples by Using ZnFe₂O₄-Based Sensing Platform: Correlation between Analyte Molecular Structure and Electronic Property of Electrode Materials, *J. Electrochem. Soc.*, 2022, **169**, 106517.
2. N. T. Anh, N. N. Huyen, N. X. Dinh, N. T. Vinh, N. Van Quy, V. D. Lam and A.-T. Le, ZnO/ZnFe₂O₄ nanocomposite-based electrochemical nanosensors for the detection of furazolidone in pork and shrimp samples: exploring the role of crystallinity, phase ratio, and heterojunction formation, *New J. Chem.*, 2022, **46**, 7090-7102.
3. D. T. N. Nga, N. L. N. Trang, V.-T. Hoang, X.-D. Ngo, P. T. Nhung, D. Q. Tri, N. D. Cuong, P. A. Tuan, T. Q. Huy and A.-T. Le, Elucidating the roles of oxygen functional groups and defect density of electrochemically exfoliated GO on the kinetic parameters towards furazolidone detection, *RSC Adv.*, 2022, **12**, 27855-27867.
4. N. T. Anh, N. X. Dinh, N. N. Huyen, P. T. L. Huong, V. N. Phan, P. D. Thang, H. Van Tuan, T. Van Tan and A.-T. Le, Promoting electron transfer kinetics and adsorption capacity for the detection of furazolidone in real food samples by using ag-core@Fe₃O₄-shell-based electrochemical sensing platform, *J. Electrochem. Soc.*, 2023, **170**, 017510.
5. A. A. Ensafi, N. Zandi-Atashbar, M. Gorgabi-Khorzoughi and B. Rezaei, Nickel-ferrite oxide decorated on reduced graphene oxide, an efficient and selective electrochemical sensor for detection of furazolidone, *IEEE Sens. J.*, 2019, **19**, 5396-5403.