Supporting Information

"Functionalized Hybrid Nanoparticles for Mercury Scavenging"

Xiaoe Li,*^a Javier Pérez-Hernández,^b Saif A. Haque,^a James R. Durrant^a and Emilio Palomares*^b

^a Department of Chemistry, Imperial College London, London, UK SW7 2AY

^b Institute of Chemical Research of Catalonia (ICIQ), Avda. Països Catalans, 16, Tarragona 43007,

Tarragona (Spain). Fax: +34 977 920 224; Tel.: +34 977 920 241; E-mail: epalomares@iciq.es.

Determination of the scavenged mercury by the mesoporous N719/TiO₂ film

The mercury scavenging experiments were performed dipping the N719/TiO₂ film (1 cm² area) in different mercury aqueous solutions ranging from 1 to 30 μ M concentrations during 10 minutes. In order to quantify colorimetrically the mercury remaining in the solution, after the scavenging process, we added 12 μ l from a 2 mM N719 stock solution in ethanol to 1 ml of the mercury aqueous solution. The UV-Visible spectrum of the mixture was measured. The results are shown in Figure S1.



Figure S1. UV-Visible spectra upon addition of 12µl from a 2 mM N719 solution in ethanol to 1 ml of HgCl₂ aqueous solution. The inset shows the mercury concentration range used.

Supplementary Material (ESI) for Journal of Materials Chemistry

This journal is (c) The Royal Society of Chemistry 2007

Using the data from **Figure S1** we can determine the mercury concentration left after the immersion of the scavenging film. To this end, we fitted the experimental data to a equation that relates the absorption changes upon mercury binding to the N719 molecule sin solution.

$\mathbf{R} = \mathbf{1.76} + \mathbf{0.03^{*}[Hg^{2+}]}$

Where $R = \frac{Abs_{\lambda 1}}{Abs_{\lambda 2}}$, $\lambda_I = 508$ nm (maximum absorbance in the absence of Hg²⁺ ions) and

 $\lambda_2 = 438$ nm (maximum absorbance in the presence of Hg²⁺ ions). The linear fitting is illustrated in **Figure S2**. Using this calibration curve we can extrapolate the mercury concentration present in the aqueous solutions after the immersion of the scavenging TiO₂/N719 films.



Figure S2. Calibration curve for the N719 solution chemosensor in ethanol.

Supplementary Material (ESI) for Journal of Materials Chemistry This journal is (c) The Royal Society of Chemistry 2007

The selectivity of the solution based mercury detection.

The colorimetric mercury sensing experiments, based upon the addition of 12μ L of N719 in ethanol to the water based samples after immersing the scavenging film, were carried out for different samples containing other metal ions as shown in **Figure S3**. The aim was to illustrate the low response of other metal ions such as Cadmium or Lead to the N719 molecules.



Figure S3. Changes on the absorption ratio (λ_1/λ_2) for several metal ions after the dipping of the TiO₂/N719 films into the aqueous solution containing different metal ions.