

Supporting information:

Synthesis and ion transport activity of oligoesters containing an environment-sensitive fluorophore

*Joanne M. Moszynski and Thomas M. Fyles**

Department of Chemistry, University of Victoria. Victoria, BC. Canada. V8W 3P6

tmf@uvic.ca

This file contains summary information on the bilayer-clamp conductance experiments conducted with the compound HO₂C-**Trip-G(E3)**-OH. The compound is abbreviated [tGE3] throughout this file.

The information consists of:

- Experimental details and comments on the structure and format of the summary files
- Five summaries of experiments conducted.

Bookmarks are inserted for the five experimental summaries.

Filename	electrolyte
JMM-3	1M CsCl
JMM-4	1M CsCl
JMM-1	1M KCl
JMM-2	1M NMe ₄ Cl
JMM-5	1M NMe ₄ Cl

Bilayer clamp experimental details.

The methodology was derived from that previously described (1). A model BC-525A bilayer clamp (Warner Instrument Corp.) was used for all planar bilayer experiments, ClampEx 8 and ClampFit 10 (Axon Instruments) were the software used for acquisition and analysis, respectively. Cups used were made of polystyrene and had 250µm diameter apertures. The lipid used in all cases was diphytanoyl phosphatidylcholine (diPhyPC) (Avanti Polar lipids). A stock solution of 25mg/mL lipid in CHCl₃ was dried under N₂ and then re-suspended in 200µL decane. For compounds that had to be pre-loaded into the lipid, 0.1 – 1mol% compound in CHCl₃ was added to the lipid mix and then dried down. The electrolytes used were 1M KCl, NaCl, CsCl, or NMe₄Cl in 10mM HEPES, 10mM TRIS, pH 7, (unadjusted).

The aperture was primed with 0.5-1µL of decane/lipid, excess solvent was removed by blowing N₂ over the aperture. The cup was then placed into the electrolyte-filled holding cell, consisting of 5mL and 3mL chambers, and salt bridges (KNO₃/Agar) and electrodes (Ag/AgCl) were attached. Bilayers were formed by brushing on 1- 1.5µL of the decane/lipid mix over the aperture, and were monitored for stability, capacitance and resistance for at least 20 minutes before test compound was added. Test compounds were added by injection from an organic solution (typically no more than 1-10µL of solution) or by breaking the lipid-only bilayer and brushing on the compound-preloaded in a lipid mix. All data were hardware filtered (8- pole Bessel filter, 1 kHz) and data was collected in a survey mode using the Gap-free protocol. Bilayers were tested repeatedly for capacitance and resistance. Once formed, 'activity' from pristine bilayers was never observed.

Conductance data summary format.




Each experiment conducted is summarized in a separate file within this portfolio. The structure of each summary is similar consisting of:

- A title of the form: file identifier {date} [compound] details of the experiment
- A plot of the applied potential as a function of the overall experiment time. Color overlays on this plot show the durations of the recorded traces, numbered sequentially 0000 – 00mn. The plot also shows when bilayer breakage occurred (if it did).
- A bar graphic showing the composition of the experiment as it proceeded. The lines establish that the indicated component was continuously present.
- Activity grid summaries for the individual traces recorded. These grids use the notation presented in ref (2). The horizontal position gives event duration in log (duration) from 1 to 10⁵ milliseconds. The vertical position gives conductance in log (conductance) from 1 to 3000 pS in 0.5 log unit steps (1,3,10,30 etc pS). The colour codes are : green = regular "square-top" behaviour; yellow = "flicker" behaviour; blue = "multiple open states" behaviour; red = "spike" behaviour; purple = "erratic" behaviour; grey indicates the region cannot be observed due to the instrument and data processing settings of the experiment .
- For each trace recorded there are one or more pages showing the overall trace and selected sections expanded to highlight the various conductance behaviours noted on the activity grid for the experiment.





(1) Eggers, P.K.; Fyles, T.M.; Mitchell, K.D.D.; Sutherland, T., *J. Org. Chem.* **2003**, *68*, 1050-1058.

(2) Chui, J.K.W.; Fyles, T.M., *Chem. Soc. Rev.* **2011** DOI: 10.1039/c1cs15099e.



JMM-3 {18-06-10} [tGE3] soln MeOH inject diPhyPC 1M CsCl 250uM styrene 'new' cup 1st try







0000





0000 continues...







0001





0001 continues...






0002







0003






JMM-4 {25-06-10} [tGE3] soln MeOH inject diPhyPC 1M CsCl 250uM styrene 'new' cup 2nd try





0000






0000 continues..






0000 continues..





0001






0003



JMM-01 {06-07-10} [tGE3] soln MeOH inject diPhyPC 1M KCl 250uM styrene 'new' cup 2nd try



broken bilayer

broken bilayer

broken bilayer


broken bilayer

broken bilayer


lipid

electrolyte


contact injection






0000






A







0002






0003




0004




0005







0006



0007




0008




0009 blank record; no data




0010




0011






0012






0013




JMM-2 {06-07-10} [tGE3] soln MeOH inject diPhyPC 1M NMe4Cl 250uM styrene 'new' cup 2nd try



0014




0015




0016 - blank file; no data

0017 - blank file; no data


0018




A




0019



0020




0021













0021 blank file; no data




0022






JMM-5 {28-06-10} [tGE3] soln MeOH inject diPhyPC 1M NMe4Cl 250uM styrene 'new' cup 1st try






0000





0000 continues..






0001





0002




0003






0004



0005





0007

