Electronic Supplementary Material (ESI) for Lab on a Chip. This journal is © The Royal Society of Chemistry 2021

## **Supplementary Information**

## An automated and portable antimicrobial susceptibility testing system for urinary tract infections

Kuo-Wei Hsu<sup>a</sup>, Wen-Bin Lee<sup>a</sup>, Huey-Ling You<sup>d</sup>, Mel S. Lee<sup>\*e</sup>, and Gwo-Bin Lee<sup>\*abc</sup>

<sup>&</sup>lt;sup>a</sup>Department of Power Mechanical Engineering, National Tsing Hua University, 30013 Hsinchu, Taiwan

<sup>&</sup>lt;sup>b</sup>Institute of Biomedical Engineering, National Tsing Hua University, Hsinchu 30013, Taiwan

<sup>&</sup>lt;sup>c</sup>Institute of NanoEngineering and Microsystems, National Tsing Hua University, Hsinchu 30013, Taiwan

<sup>&</sup>lt;sup>d</sup>Laboratory Medicine, Kaohsiung Chang Gung Memorial Hospital, Chang Gung University, Kaohsiung 83301, Taiwan

<sup>&</sup>lt;sup>e</sup>Department of Orthopaedic Surgery, Kaohsiung Chang Gung Memorial Hospital, Chang Gung University, Kaohsiung 83301, Taiwan.

Table S1: Comparison of Fluidigm® C1 system and our devices

	Fluidigm® C1 system	This study
Microvalve	Normally-open	Normally-closed
Micropump	Peristaltic pumping	One-time injection
Micromixer	None	Membrane-type