

In situ Raman monitoring of electroreductive dehalogenation of aryl halide at Ag/ aqueous solution interface

Ying-Hua Zhou ^{a,b,#}, Chen-Chen Jiang ^{b,#}, Zhou Yu^{b,#}, Ya-Hao Wang ^{b,*}, Ju-Fang Zheng ^{b,*}, Xiao-Shun
Zhou ^{b,*}

a. Jinhua Education College, Jinhua 321000, China

b. Key Laboratory of the Ministry of Education for Advanced Catalysis Materials, Institute of
Physical Chemistry, Zhejiang Normal University, Jinhua 321004.

#. Corresponding authors

E-mails: jfzheng@zjnu.cn; yahaowang@zjnu.edu.cn; xszhou@zjnu.edu.cn

Supplementary Figures




Fig. S1 SHINERS spectra of commercial TP and BPDT on the Ag electrode compared to CTP and BTP at -0.4 V.