




SEI_I: TG/DTG trends registered for the samples reacted at 30, 45, 60, and 70°C in 10M NaOH solutions.



SEI_II: TG/DTG trends registered for the samples reacted at 30, 45, 60, and 70°C in 12M NaOH solutions.



SEI_III: TG/DTG trends registered for the samples reacted at 30, 45, 60, and 70°C in 15M NaOH solutions.



SEI_IV: TG/DTG trends registered for the samples reacted at 30, 45, 60, and 70°C in 17M NaOH solutions.